
MPV TRUMA

Calimero® plus

SIDS-Monitor

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Special notes

1. The monitor is not defibrillator-proof. Disconnect the monitor from the patient before using a defibrillator.
2. The monitor is not designated for use in explosion-endangered locations.
3. Do not use the monitor simultaneously with other monitors which monitor respiration by impedance measurement (measuring current passing through the patient's body). This may damage the Calimero®plus SIDS monitor and / or disrupt the monitoring process.

Responsibility of the manufacturer

The manufacturer and importer are responsible for the effects on safety, reliability and performance only when:

1. assembly / installation, extensions, readjustments, modifications or repairs are carried out only by personnel authorized by the manufacturer.
2. Calimero®plus SIDS monitor is used in accordance with the instruction manual.
3. the electrical installation at the location complies with VDE / IEC regulations.

No liability is accepted for repairs, modifications, adjustments or similar interventions carried out by the user's technical personnel.

This also applies to devices which the manufacturer has designated as "reparable by the customer" and for which appropriate documents are available.

Brief instructions

Please read the instructions for use carefully before first use of your Calimero®plus SIDS-Monitor. The following brief instructions will help you to find the most important control elements quickly.

Operating

1. Attach the disposable electrodes to your child's chest as in the diagram (see Fig. 1), ensuring that the skin is dry and oil-free.
2. Assign the electrodes to the colors of the patient main cable. Red and yellow may be exchanged, the blue body sensor must be positioned at the location indicated (see also Fig. 1).
Ensure that the body sensor is not kinked!
3. Use the enclosed Leukosilk or Leukopor adhesive tape to attach the body sensor (see Fig. 1). One strip across the sensor is enough. Connect the body sensor to the protruding black socket on the ECG cable.
After use, remove the adhesive tape with great care. Carefully remove any adhesive material left on the sensor.
4. Connect the patient main cable to the socket ⑩ "Socket for patient main cable".
5. Press the ON button ⑧ to start the SIDS-Monitor and keep it pressed for about 1 second.

Diagram of electrode connections


Fig. 1


Fig. 2

Other possible electrode positions, e.g. in the event of frequent false alarms (bradycardia).

Purpose and function

Dear Customer,

With Calimero®plus SIDS monitor you have purchased an easy-to-use device for monitoring both the respiratory rate and heart rate of your child.

Calimero®plus SIDS monitor monitors both the respiratory rate and heart rate of your child while he is sleeping, and gives an alarm in the event of a dangerous interruption to respiration, or if the heart rate exceeds maximum or minimum limits.

Calimero®plus SIDS monitor saves all processes automatically to an accuracy of seconds (ON / OFF, RESET, all types of alarms, heart and respiratory rate before and after an alarm). No special operating procedures are required for this function.

Before you start using the SIDS monitor, please note the following:

All persons who take care of your child should carefully read through the instructions for use. This will prevent operating errors and misunderstandings.

These instructions for use represent an integral part of the device and should always be available wherever it is used.

With Calimero®plus SIDS monitor, you only need three electrodes – one for respiration and two for cardiac monitoring. An additional electrode may be needed if the heart monitoring malfunctions due to external electromagnetic influences.

For special cases, it is also possible to use the device only for respiration or only for heart rate monitoring. Consult your doctor and contact the retailer or manufacturer concerning the operating procedures.

In the case of an alarm

Make a brief check to determine the nature of the alarm:

- Respiration alarm? (Red LED ④ lights up)
- Heart rate alarm?
(One or both red LEDs ② + ③ light up)
(bradycardia, tachycardia).
Switch off the SIDS monitor and proceed as described in the section “What should you do in an alarm situation?” on Page 37.
- Technical alarm: red LED ⑨, ⑩ or ⑪ lights up
(battery, sensor for respiration or heart).
Switch off and proceed as described in “Notes on faults” on Page 37.

Operating instructions

Installing the device

Place the device in a position where the displays are easily visible on entering the room. In addition, the equipment should not be placed directly above the bed. Otherwise, the device could fall down on the bed if the child pulls the cable.

Before switching on the Calimero[®]plus SIDS monitor, insert the plug of the patient main cable in the socket ⑥ at the top of the device.

It is best to bring out the cable at the bottom of the bed to prevent your child from using it as a toy. After positioning the sensors (see Page 29), press button ③ to switch on the device.

As soon as it is switched on, the device performs a self-test for all function displays: all monitoring displays light up in turn. A short self-test tone also sounds. This allows you to test the correct functioning of all displays each time the equipment is switched on. The flashing green monitoring lamp ⑤ together with a soft clicking tone on the Calimero[®]plus SIDS monitor indicate regular respiration. You can enable / disable the monitoring tone by pressing the "Respiration monitoring tone" button ⑥. Regular heart activity is monitored by the "Heart function display" ⑫.

Memory function

The memory function is always enabled. No special operating procedures are required for this function.

Alarm delay

Babies occasionally experience interruptions in respiration of more than 10 seconds. This is not unusual.

These irregularities in the respiratory rhythm decline with increasing age and, in most cases, they have disappeared by three years of age.

If, however, the respiratory arrest exceeds the above specified period and /or occurs more frequently, then your child may be at risk.

To compensate for individual fluctuations, the delay time after which the alarm is given can be adjusted with the rotary switch ③.

First set a period of 15 seconds (for babies) or 10 sec. (from approximately 2 months of age).

If the alarm occurs frequently although there is no reason for concern, the alarm delay time can be increased after consulting the attending pediatrician. (20 or 25 seconds)

Respiration alarm

In an alarm situation, the red control lamp ④ lights up and a loud intermittent signal tone sounds.

The alarm is triggered in the absence of an impulse from the sensor = emergency situation = interrupted respiration. In the event of an alarm situation, refer to the information given under "What should you do in an alarm situation?" on Page 40.

Setting limits / Testing

Cardiac monitoring - setting the limits

The rotary switches ① and ⑭ set the alarm limits in heart rate monitoring. Since these limits vary according to the age of the child, you should consult your pediatrician to discuss these settings. For recommended values, see table below:

Age	Heartrate	possible settings	
		min	max
Premature infant	120 - 170	90	240
Newborn	115 - 150	90	220
6 months	100 - 140	80	220
1 year	100 - 140	70	200
2 years	80 - 130	60	200

When adjusting min./max. values, turn the corresponding switch until it latches. Do not set intermediate values.

A decrease in the heart rate (bradycardia) of a baby is much more serious than an increase in heart rate (tachycardia), since bradycardia is often caused by a lack of oxygen, which is a life-threatening situation and needs immediate treatment.

Testing

Switch on the device once a week with the cables and sensors connected, but without attaching the electrodes to your child.

The SIDS monitor must respond with a definite alarm.

Battery alarm / battery exchange / Cleaning

Battery alarm / battery exchange

The Calimero®plus SIDS monitor plus is operated with a 9 V battery (alkaline battery if possible) and with average use can operate on this battery for about 1,5 months. This operating time can be extended with a special lithium battery to a period of between 2 and 3 months. Please contact the manufacturer or your dealer in this respect.

Before the actual battery alarm, the SIDS monitor gives a battery warning:

The yellow LED ⑤ lights up every 10 sec. during operation. In this case, the battery must be replaced within the next 1 to 2 days. You can suppress the flashing for a limited period of time by pressing the RESET button ⑦. If the battery is further discharged during operation, the actual battery alarm will be activated. A continuous warning tone sounds and the yellow "Battery alarm" indicator ⑥ lights up constantly.

In this case, the battery must be replaced immediately.

Please note that acid leaking from a discharged battery could cause damage to the device. To exchange the battery, open the battery compartment ⑧ and insert the new battery as shown in the drawing.

Advice on the disposal of batteries:

After use, batteries (in particular alkaline batteries and nickel cadmium batteries) must be disposed of as hazardous waste or at your dealer's collection center. Do not include it in domestic waste!

Battery alarm / battery exchange / Cleaning

Cleaning

If required, carefully wipe the case with a slightly damp cloth. Make sure that no moisture enters the device.

Malfunctions

- Problem:** The device gives an alarm immediately after it is switched on. (10 oder 11)
- Remedy:** Check sensor plug and sensors.
- Problem:** The device gives an alarm and the yellow lamp 9 lights up.
- Remedy:** Fit new battery.
- Problem:** The device gives a warning tone and the red LED 11 lights up.
- Remedy:** Check electrodes for proper contact. If necessary fit new disposable electrodes. Apply electrodes to a new position (Fig. 2).
- Problem:** The device gives a warning tone and the red LED 10 lights up.
- Remedy:** Check the connection or exchange the body sensor. Frequent unexplained apnoe alarms, particularly in babies aged 1 – 4 weeks.
- Remedy:** Set the slide switch in the battery compartment to "E" (sensitive). After breathing is stabilized, reset to "N" (normal)! In babies from 6 weeks old, ensure that your monitor is in the "N" position.

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Technical data

Accessory of other manufacturers (disposable electrodes) approved for the Calimero®plus SIDS monitor:

Manufacturer	Description
Kendall	H87V (disposable electrodes)
Medicotest A/S.DK	NF-60-K/HC (disposable electrodes)

Scope of delivery:	Calimero®plus SIDS monitor, body sensor, 12 disposable electrodes, patient main cable, 9 V battery, instructions for use, training video, case.
Operation mode:	Battery-operated, 9 V monobloc battery, Battery type: alkaline or lithium
Battery replacement:	after approx. 1,5 or 3 months if operated with lithium battery
Operation temperature:	10 – 40° C
Storage temperature:	4 – 50° C
Function displays:	Respiration signal: green LED 5 Heart signal: green LED 12 Respiration monitoring tone: interruptible
Alarms:	Interrupted respiration: red LED 4 and warning tone Bradycardia: red LED 2 and warning tone Tachycardia: red LED 13 and warning tone Discharged battery: yellow LED 9 and tone Defective cable or electrodes incorrectly applied: red LEDs 10 and 11 and tone

Technical data

Switches:	Rotary switch for min. and max. heart rates, Rotary switch for respiration alarm delay, Switch for respiration sensitivity (in battery compartment)
Alarm delay times:	10, 15, 20, 25 sec. (respiration)
Min. / max. values:	60, 70, 80, 90 / 180, 200, 220, 240 BPM (heart rate)
Connections:	Patient main cable, external signal transmitter (optionally hospital pager or alarm center)
Dimensions:	l x w x h (in mm): 152 x 78 x 29
Weight:	Approx. 245 g
Protection class:	BF
Classification to MPG:	IIB, 93 / 42 / EWG
 0125:	Calimero® plus SIDS monitor is conform to the requirements of Guideline 93/42/EWG and Medical Products Act
Manufacturer:	Schulte Elektronik, Olsberg

Technical details subject to change.

What should you do in an alarm situation?

Important note

Before using this equipment, ensure that you make yourself familiar with resuscitation measures. Without this knowledge you cannot provide help in the event of an emergency. To make yourself familiar with the necessary procedures, we recommend you to read the relevant literature in addition to taking part in practical training under the direction of a trained person. For further training possibilities ask your physician or contact your local public health authority.

The following information represents only a very short summary of the measures which must be taken in the event of an emergency:

- Switch off the device and check the following points:
 - Has your child stopped breathing?
 - Has his skin turned blue? (lips, face, fingers)
 - Is his pulse slow or is there no pulse at all?
- If one of the above mentioned points is applies, ask another person to call for the emergency physician immediately.
- Always start with resuscitation measures immediately.
- Take your child out of bed and try to wake him up by shaking him.
- If your child does not start breathing spontaneously, start giving artificial respiration immediately. If you cannot feel any pulse, start with cardiac massage immediately.
- Continue with these measures until the emergency physician arrives.

We must point out once more that the action under point 5 can only be taught by a trained person. In addition, knowledge of resuscitation measures for adults is not sufficient since there are a number specific differences involved. The video film "A chance for life" contains important information on resuscitation.

Literature

Video film (VHS): „wohlbehütet“ (A chance for life“), first-aid measures for babies and infants.

Installation instructions for the read-out system

Installing the programm


You will find a detailed description in the program under "HELP"

- Switch off your Computer
- Connect the PADMOUSE® system to an unassigned serial port (COM 1-4) of your computer.
- Start your Computer
- Insert installation disc in 3,5" disk drive.
- Select <setup.exe> on disk.
- Follow the instructions of the installation program.
- Start the <SISS> monitor evaluation software by clicking on the program icon, or select from the start menu.
- For data transfer follow the instructions on the display.
- The PADMOUSE® system can be permanently connected to the computer.

System requirements:

- PC with operating system Windows 3.1x/ 95/ NT
- 1 unassigned serial port (COM 1-4)
- 2 MB available main memory
- Pad of the PADMOUSE® system

PC with windows
3.1x/ 95/ ME/
NT/ 2000

Connection to
serial port
(COM 1-4)

Data read-out

1. Ausleseprogramm - Switch on PC and double-click onto "SISS" to start the read-out program;
- Click once to select "Data transfer".
2. Uhrzeit und Datum Check the time on your PC and correct if necessary. The correct time and date is essential to obtain correct data; confirm time / date.
3. Monitor aufsetzen: The system dispenses with the need for a cable connection between the device and PC. Simply position Calimero®plus SIDS monitor on the pad connected to the computer.
4. Auslesen starten: Press " " and "RESET" simultaneously. During transfer, the red LED "ALARM" flashes after ... SECONDS". The transfer process takes 5 - 60 sec. depending on the number of events.
5. Auslesen beenden: The PC informs you as soon as the transfer is completed.
Remove the device from the pad.
6. Daten im Monitor After transfer you can clear the memory in the device by pressing RESET. The SIDS monitor switches off automatically. If the data are not deleted, the device switches off automatically after 2 minutes.
7. Daten auswerten: Click on an event. The screen shows additional information relating to the event selected. Click on "Evaluation" to view the statistics of events. Select "Print" to print out the values displayed.

Information on battery replacement with the Calimero®plus SIDS monitor: The time limit for removing a discharged battery and inserting a new battery is 5 minutes. If you do not comply with this limit, stored data could be lost.

Windows and Windows 95/98/ME/NT/2000 are registered trade marks of Microsoft Corporation.

Order Details

Description	Part No.	VE
Calimero®plus SIDS monitor for breathing and heart monitoring with data store, complete with body sensor for breathing, single use electrodes for heart monitoring (packet of 12), patient main cable, instructional video in carry case and instruction manual	M 72600-00	1
Accessories – Data analysis:		
Calimero®plus analysis system complete with software for data analysis and editing as well as Pad-Mouse for data transfer to PC	M 70060-02000	1
Calimero®plus instructional video "Risk factors and recommendations in relation to sudden infant death syndrome with all first aid measures for babies and infants"	M 70060-02500	1

Order Details

Description	Part No.	VE
Accessories – Replacement parts:		
External signaller for Calimero®plus SIDS-Monitor with 15 m connection cable, 9 V battery	M 70060-03000	1
Patient main cable	M 70060-03300	1
Calimero®plus body sensor for breathing	M 70060-00400	1
Electrode for heart monitoring (packet of 30)	M 70060-04200	30
Long-life electrode belt for heart monitoring of children with 2 longlife electrodes	M 70060-04300	1
Battery storage lid Calimero®plus	M 70060-04500	1
Calimero®plus carry case	M 70060-03900	1
Calimero®plus instruction manual	M 70060-00500	1
Maintenance – Update:		
Calimero®plus device servicing Upgrade to latest technical level, software update and functional test	M 70060-01100	1
Accessories – First Aid:		
Resuscitator mask-kit I for infants up to 10 kg body weight, made of silicone with silicone masks size 0 and size 1 and patient valve	M 84800-00	1
Resuscitator mask size 0 for infant small, 100% silicon	M 80070-01000	1
Resuscitator mask size 1 for infant large, 100% silicon	M 80070-01100	1

Warranty

This device had been manufactured with due care and had been thoroughly checked before it left our factory.

We therefore provide a warranty for this product for a period of 1 year from the date of purchase according to the following conditions:

Within the warranty period, we will remedy any deficiencies in the device free of charge which are caused by material or manufacturing defects, at our option either by repair or replacement of parts or of the device. Replaced parts and devices pass into our ownership. This warranty does not cover damage to or deficiencies in the device resulting from improper use (e.g. operation with the wrong current or voltage, connection to unsuitable power sources, breakage), or if the defect has only an insignificant effect on the value or utilization of the device.

The manufacturer shall be responsible only for the effects on safety, reliability and performance of the equipment if:

- expansions, readjustments, modifications or repairs have been carried out by personnel authorized by the manufacturer and
- the electrical installation at the location conforms to the requirements specified by IEC regulations and
- the device is used in accordance with the instruction manual.

Remedial action taken under the warranty provisions does not affect the warranty period. In the event of claims under warranty, return the complete device with any accessories correctly packed along with your receipt to the manufacturer whose address is specified below, or to your dealer. All goods must be returned carriage-free. If postage is unpaid or insufficient, we will refuse to accept the delivery.

This warranty is valid only if the date of purchase is confirmed by the dealer's stamp and signature. If you bought the device directly from the manufacturer, you need only indicate the invoice number.

Any other claims, in particular claims for damages or compensation, are excluded. This also applies to those cases where the defect has obviously been caused by a technical failure of the device.

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