AIR FORCE ONE

CE 0123

DE AIR FORCE ONE KOMPRESSOR-VERNEBLER

EN AIR FORCE ONE AEROSOL THERAPY SYSTEM

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






G

H

- | | | |
|-----------------------|------------------------------|----------------------|
| 1. Compressor housing | 7. Air-vent openings | 13. Mask (PP+TPE) |
| 2. ON/OFF switch | 8. Handle | 14. Spare air ers |
| 3. Power cord | 9. Air tube (PVC) | 15. Accessories case |
| 4. Air outlet | 10. Nebulizer (PP+Silicone) | 16. Carrying case |
| 5. Filter case | 11. Mouthpiece (PP+Silicone) | |
| 6. Nebulizer holder | 12. Nosepiece (PP) | |

List of Material

Part	Material
Nebuliser	Polypropylene
Mouthpiece	Polypropylene
Adultmask and Childmask	Polypropylene and Thermoplastic

Read carefully the Instructions for use before Operation.

WARNINGS

1. Use the device only as described in this manual and therefore as an aerosoltherapy system, following the indications of your doctor. Any use different from the intended one is to be considered improper and hence dangerous; the manufacturer cannot be held liable for any damage caused by improper, incorrect and/or unreasonable use, or if the equipment is connected to electrical installations which do not comply with current safety regulations.
2. Keep this manual for future reference.
3. Do not operate the unit in presence of any anaesthetic mixture in ammable with oxygen or nitrogen protoxide.
4. The correct functioning of the equipment can be a ected by electromagnetic interferences which exceed the limits indicated by the European standards in force. In case this device interferes with other electrical devices, move it and plug it to a di erent power socket.
5. In case of failure and/or malfunction, read the "POSSIBLE PROBLEMS AND HOW TO SOLVE THEM" section. Do not handle or open the compressor housing.
6. For repair operations address only to a technical service centre authorized by the manufacturer and require the use of original spare parts. The nonobservation of the above mentioned indications can compromise the device safety.

7. Comply with the safety regulations concerning the electrical devices and in particular:
 - use only original accessories and components;
 - never submerge the unit in water;
 - never wet the device, it is not protected against water penetration;
 - never touch the unit with wet or moist hands;
 - do not leave the unit exposed to the weather elements;
 - place the unit on a stable and horizontal surface during its operation;
 - the use of this device by children and disabled requires always the close supervision by an adult with full mental faculties;
 - do not pull the power cord or the device itself to unplug it from the power socket;
 - the power plug is the separation element from the grid power; keep the plug accessible when the device is in use.
8. Before plugging in the device, make sure that the electrical rating, shown on the rating plate on the bottom of the unit, corresponds to the mains rating.
9. In case the power plug provided with the device does not fit your wall socket, address to qualified personnel for the replacement of the plug with a suitable one. In general, the use of adapters, simple or multiple, and/or extension cables is not recommended. If their use is indispensable, it is necessary to use types complying with safety regulations, paying however attention that they do not exceed the maximum power limits, indicated on adapters and extension cables.
10. Do not leave the unit plugged in when not in use; unplug the device from the wall socket when it is not operated.
11. The installation must be carried out according to the instructions of the manufacturer. An improper installation can cause damage to persons, animals or things, for which the manufacturer cannot be held responsible.
12. The power cord of this device cannot be replaced by the user. In case of a power cord damage, address to a technical service centre authorized by the manufacturer for its replacement.
13. The power supply cord should always be fully unwound in order to prevent dangerous overheating.
14. Before performing any maintenance or cleaning operation, turn off the device and disconnect the plug from the main supply.
15. Some parts of the unit are so small that they may be swallowed by children; keep the equipment out from children's reach.
16. If you decide not to use the device any longer, it is recommended to dispose of it according to the current regulations.
17. Make sure to:
 - use this device only with medicines prescribed by your doctor;
 - make the treatment using only the accessory recommended by your doctor depending on the pathology ;
 - use the nosepiece accessory only if expressly indicated by your doctor and paying attention to NEVER introduce the bifurcations in the nose, but only bringing them as close as possible.
18. Check in the medicine package leaflet for possible contraindications for use with common aerosoltherapy systems.
19. To avoid strangulation and entanglement, keep cable and air tubes out of reach of young children.
20. Do not position the equipment so that it is difficult to operate the disconnection device.

21. For greater hygienic safety, we recommend you to avoid using the same accessories for more than one person.

Instructions for use

The device must be checked before each use, in order to detect possible functioning anomalies and/or damages due to transport and/or storage.

During inhalation, sit upright and relaxed at a table and not in an armchair, in order to avoid compressing your respiratory airways and impairing the treatment effectiveness.

1. After unpacking the device, check it for visible damages or defects; pay particular attention to cracks in the plastic housing, which may expose electrical components. Check for accessories integrity.
2. Before using the device, proceed with the cleaning operations as described in the "CLEANING AND MAINTENANCE" section.
3. Open the nebulizer by turning counterclockwise the top (Pic. B).
4. Make sure that the medicine conduction cone is properly fitted on the air conduction cone inside the nebulizer (Pic. C).
5. Put the prescribed quantity of medicine into the nebulizer (Pic. D).
6. Close the nebulizer by turning clockwise the two parts, paying attention that they are thoroughly sealed (Pic. E).
7. Connect one end of the air tube to the nebulizer (Pic. F) and the other end to the air outlet (Pic. G) on the device.
8. If you need to use the mask, insert it directly onto the nebulizer (Pic. H).
9. Plug the device into the wall socket, making sure that the main supply corresponds to the electrical rating of the device.
10. To start the treatment, set the I/O switch into the "I" position.

ATTENTION: this device is for intermittent use: 30 min ON / 30 min OFF.

11. Inhale the aerosol solution using the prescribed accessory.
12. When the treatment has been completed, switch the unit off by setting the I/O into the "O" position and disconnect the plug from the wall socket.
13. Wash the nebulizer and its accessories as described in the "CLEANING AND DISINFECTION OF THE ACCESSORIES" section.
14. The device requires no calibration.

Cleaning and maintenance

The cleaning of the device must be carried out by using a soft and dry cloth and non-abrasive cleansers.

ATTENTION: during cleaning operations, make sure the internal parts of the unit are not in contact with liquids and that the power plug is disconnected.

Cleaning and disinfection of the Accessories

Follow carefully the cleaning and disinfection instructions of the accessories as they are very important for the device performances and the therapy success.

After each treatment:

Disassemble the nebulizer by turning counterclockwise the top and remove the medicine conduction cone. Wash the components of the disassembled nebulizer, the mouthpiece, nosepiece and mask by using tap water; dip in boiling water for 5 minutes. Reassemble the nebulizer components and connect it to the air-outlet (Pic.G), switch the device on and let it work for 10-15 minutes. The air tube must be washed with warm water. Do not boil nor autoclave the air tube.

Replacement of the nebulizer

The nebulizer must be replaced after a long period of inactivity, in case it shows deformations or breakings, or when the nebulizer nozzle is obstructed by dry medicine, dust, etc. DeVilbiss Healthcare recommends to replace the nebulizer after a period between 6 months and 1 year depending on the usage.

Use original nebulizers only.

Replacement of the air filter

In normal usage conditions, the air filter must be replaced approximately after 500 working hours or after each year. DeVilbiss Healthcare recommends to periodically check the air filter (10 - 12 treatments) and, if the filter shows a grey or brown colour or is wet, replace it. Extract the filter (Pic. A-5) and replace it with a new one. Do not try to clean the filter for reusing it. The air filter shall not be serviced or maintained while in use with a patient.

Use original filters only. Do not use the device without filter.

TROUBLESHOOTING

Contact an authorized Customer Service Centre in case you need help about the use and the maintenance of the equipment.

PROBLEMS:

1. The device does not switch on
2. The device does not nebulize or nebulizes weakly

SOLUTIONS:

1. Make sure the power plug is firmly fitted to the wall socket;
2. Make sure that the device has been operating within operating limits indicated in this manual (30 min ON / 30 min OFF).
3. Make sure that the ends of the air tube are tightly fitted to the main unit and nebulizer;
4. Verify that the nebulizer is not empty or has been filled with the proper quantity of medicine (max 8 ml). Verify that the nebulizer nozzle is not obstructed.

Should the device not start working properly again, address to **DEVILBISS CUSTOMER SERVICE**

MAINTENANCE AND REPAIRS

In case of failure, address to qualified personnel authorized by DeVilbiss Healthcare Do not open the device in any case. The unit has no user-serviceable parts within and does not need internal maintenance or lubrication.

TECHNICAL SPECIFICATIONS

Power supply:	230V 50 Hz AC
Minimum filling volume:	2 ml
Maximum filling volume:	8 ml
Nebulization rate:	0,3 ml/min
Maximum pressure:	2,2 bar
Weight:	1,5 kg
Size:	200 x 145 x 97 mm
Noise level:	52 dBA
Expected service life:	1000 hours

- Class II device as regards protection against electric shocks.
- Nebulizer, mouthpiece and masks are type BF applied parts.
- Device not protected against sprinkles.
- Device not suitable for use in presence of anaesthetic mixtures inflammable with air, oxygen or nitrogen protoxide.
- Device for intermittent use (30 min ON / 30 min OFF)
- Device not suitable for use in anaesthesia or lung ventilation systems.

Aerosol characteristics according to EN 13544-1 ANNEX CC

Flow:	5,1 l/min
Aerosol output:	0,45 ml
Aerosol output rate:	0,02 ml/min
Particle size (MMAD):	ca. 5 µm

NOTE: The technical specifications may change without prior notice.

Applied standards Electric Safety Standards:

Electric Safety Standards EN 60601-1. Electromagnetic Compatibility according to EN 60601-1-2. Class IIa Medical Device according to the 93/42/EEC European Directive on Medical Devices

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







EN WARRANTY CONDITIONS

The warranty consists in the replacement and/or free repair of originally defective components. The warranty does not cover the accessories supplied with the device and those parts subject to normal wear and tear. The device must be repaired by authorized technical service centres only. Transport costs for the device are at user's charge.

Any repair out of warranty conditions is at user's charge. The warranty is void if the device has been tampered, if the defect results from improper use or if the damage is not imputable to the manufacturer. The warranty does not involve any compensation for any kind of damages, either direct or indirect, to persons or things, occurred when the device is malfunctioning. The warranty is valid from the date of purchase certified by the purchase ticket or invoice which must be peremptorily attached to the warranty coupon.

The lack of the warranty coupon properly filled in and validated by a purchase certification voids the warranty.

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COUPON TO BE RETURNED IN CASE OF REPAIRS

EN

The warranty is valid only if accompanied by the purchase Ticket/Invoice.

PRODUCT TYPE

Model:

Serie No.:

Date of purchase:

BUYER'S DATA

Name:

Address:

Phone Number:

Defect description:

Signature for acceptance of the above mentioned warranty conditions.

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EN DEVICE SYMBOLS

I	ON		Manufacturer
O	OFF		Serial Number
~	Alternating current		Keep dry
	Class II		Attention
	Type BF	RH	Air humidity
CE 123	Complying with the 93/42/CEE Directive		Temperature
P	Atmospheric pressure		Refer to instructions

Correct Disposal of This Product (Waste Electrical and Electronic Equipment) (Applicable in the European Union countries and countries with separate collection systems) This symbol on the product or its literature indicates that the product complies with the electrical and electronic equipment and should not be disposed of with household waste. The user is responsible for assigning the end of life, to the appropriate collection, subject to the penalties provided for in current legislation on waste. For further information regarding the collection systems available contact your local waste disposal service.

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Globalcare Medical Technology Co., Ltd.
Donawa Lifescience Consulting Srl
DeVilbiss Healthcare GmbH